

LIGHT *Sheer*[™]

Diode Laser Systems

| **Clinical Results**

Permanent Hair Reduction—Proof of Results

CHIN

Fitzpatrick Skin Type IV
46 months after last treatment

Fluence

2 treatments at 30 J/cm²

2 treatments at 40 J/cm²

Pulse Duration

2 treatments at 15 ms

2 treatment at 20 ms



ARM

Fitzpatrick Skin Type II
14 months after last treatment

Fluence

40 J/cm²

Pulse Duration

20 ms

Number of Treatments

2



LEG

Fitzpatrick Skin Type III
24 months after last treatment

Fluence

40 J/cm²

Pulse Duration

20 ms

Number of Treatments

3



- FDA cleared for permanent hair reduction on all skin types I - VI including tanned skin
- Proven success in the treatment of pseudofolliculitis barbae (PFB)
- Long-term clinically documented studies on permanent hair reduction

Laser Hair Removal


FACE

Fritzpatrick Skin Type V
Four months after last treatment
Fluence: 22-24 J/cm²
Pulse Duration:
Treatments 1-4: 30 ms
Treatment 5-7: Auto



CHIN

Fritzpatrick Skin Type III
Thirty-six months after last treatment
Fluence 5 treatments at 40 J/cm²
Pulse Duration
5 treatments at 20 ms



AXILLA

Fritzpatrick Skin Type V
Four months after last treatment
Fluence: 22-24 J/cm²
Pulse Duration:
Treatments 1-3: 30 ms



Gold-Standard LightSheer™ Diode Laser Technology for Hair Removal

LIGHTSheer™


UPPER BACK and SHOULDERS

Fritzpatrick Skin Type V
Five months after last treatment
Fluence: 23 J/cm²
Pulse Duration:
Treatments 1-3: 30 ms
Treatment 4: Auto



LOWER BACK and SACRAL REGION

Fritzpatrick Skin Type V
Five months after last treatment
Fluence: 22 J/cm²
Pulse Duration:
Treatments 1-4: 30 ms



The LightSheer 100 ms and 400 ms pulse duration was not available at the time of these patient treatments.

Lumenis Global Sales Offices

Lumenis Inc.
2400 Condensa Street
Santa Clara, CA 95051 USA
Tel 408 764 3000
877 586 3647
Fax 408 764 3999
800 505 1133

Lumenis (HOLLAND) BV
Amstelveen, The Netherlands
Tel + 31 20 347 50 60
Fax + 31 20 347 50 99

Lumenis (GERMANY) GmbH
Dreieich-Dreieichenhain,
Germany
Tel + 49 6071 922 0
Fax + 49 6071 922 400

Lumenis (FRANCE) SARL
Orsay, France
Tel + 33 1 69 33 14 20
Fax + 33 1 60 19 57 23

Lumenis (FRANCE) SARL
Orsay, France
Tel + 33 1 69 33 14 20
Fax + 33 1 60 19 57 23

Lumenis (ITALY) Srl.
Rome, Italy
Tel + 39 06 90 75 230
Fax + 39 06 90 75 269

Lumenis (UK) Ltd.
London, UK
Tel + 44 20 8324 4200
Fax + 44 20 8324 4222

Lumenis (JAPAN) Co. Ltd.
Tokyo, Japan
Tel + 81 3 5789 8300
Fax + 81 3 5789 8310

Lumenis (Holland) BV
Beijing Office, Beijing, China
Tel + 86 10 6510 2620
Fax + 86 10 6510 2621

Lumenis (HONG KONG) Ltd.
Hong Kong
Tel + 852 2174 2800
Fax + 852 2722 5151

Internet: www.lumenis.com
Email: information@lumenis.com

Photo Credit:

*Christine C. Dierickx, M.D.,
Valerie B. Campos, M.D.,
R. Rox anderson, M.D.,
Wellman Laboratories of Photomedicine,
Harvard Medical School, Boston, MA
Dr. Nimish Patel, M.D. Plastic Surgeon
The Laser Center, Ahmedabad, India*



Lumenis, its logo and LightSheer are trademarks of the Lumenis Group of Companies. ©2005, the Lumenis Group of Companies. All rights reserved.